

Product datasheet for **GA103577**

Human PKNOX1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PKNOX1 CRISPRa kit - CRISPR gene activation of human PBX/knotted 1 homeobox 1
Symbol:	PKNOX1
Components:	GA103577G1, PKNOX1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103577G2, PKNOX1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103577G3, PKNOX1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P55347</u>
Synonyms:	pkonx1c; PREP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)