

## Product datasheet for **GA103592**

### Human Phospholipase C beta 4 (PLCB4) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | PLCB4 CRISPRa kit - CRISPR gene activation of human phospholipase C beta 4 |
| Symbol:       | PLCB4  |
| UniProt ID:   | <a href="#">Q15147</a>   |
| Synonyms:     | ARCND2; PI-PLC   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)    |

#### Product images:



[View online >](#)