

## Product datasheet for **GA103594**

### Human PLCL1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PLCL1 CRISPRa kit - CRISPR gene activation of human phospholipase C like 1 (inactive)
Symbol:	PLCL1
UniProt ID:	<a href="#">Q15111</a>
Synonyms:	MGC126580; MGC138190; OTTHUMP00000163640; phospholipase C, epsilon; phospholipase C-like 1; PLC-L; PLCE; PLCL; PLDL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)