

Product datasheet for **GA103597**

Human Phospholipase D1 (PLD1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLD1 CRISPRa kit - CRISPR gene activation of human phospholipase D1
Symbol:	PLD1
UniProt ID:	Q13393
Synonyms:	Phospholipase D gene 1 probably with two splicing variants A and B; phospholipase D1; phospholipase D gene 1, probably with two splicing variants A and B; Pld1a; Pld1b; Plda; Pldb
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)