

## Product datasheet for **GA103601**

### Human Pleckstrin (PLEK) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PLEK CRISPRa kit - CRISPR gene activation of human pleckstrin
Symbol:	PLEK
UniProt ID:	<a href="#">P08567</a>
Synonyms:	P47
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)