

Product datasheet for **GA103608**

Human FXYD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FXYD3 CRISPRa kit - CRISPR gene activation of human FXYD domain containing ion transport regulator 3
Symbol:	FXYD3
UniProt ID:	Q14802
Synonyms:	MAT8; PLML
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)