

Product datasheet for **GA103615**

Human PLS1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLS1 CRISPRa kit - CRISPR gene activation of human plastin 1
Symbol:	PLS1
Components:	GA103615G1, PLS1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103615G2, PLS1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103615G3, PLS1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14651</u>
Synonyms:	I-PLASTIN; I isoform; Plastin-1; plastin 1; plastin 1 (I isoform)
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)