

Product datasheet for **GA103636**

Human Pancreatic Lipase (PNLIP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PNLIP CRISPRa kit - CRISPR gene activation of human pancreatic lipase
Symbol:	PNLIP
UniProt ID:	P16233
Synonyms:	PL; PNLIPD; PTL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)