

## Product datasheet for **GA103638**

### Human Pancreatic Lipase Related Protein 2 (PNLIPRP2) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PNLIPRP2 CRISPRa kit - CRISPR gene activation of human pancreatic lipase related protein 2 (gene/pseudogene)
Symbol:	PNLIPRP2
UniProt ID:	<a href="#">P54317</a>
Synonyms:	PLRP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)