

Product datasheet for **GA103655**

Human POLR2C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	POLR2C CRISPRa kit - CRISPR gene activation of human RNA polymerase II subunit C
Symbol:	POLR2C
UniProt ID:	<u>P19387</u>
Synonyms:	DNA directed RNA polymerase II polypeptide C; hRPB33; hsRPB3; polymerase (RNA) II (DNA directed) polypeptide C (33kD); polymerase (RNA) II (DNA directed) polypeptide C, 33kDa; RNA polymerase II subunit 3; RPB3; RPB3, RPB31, hRPB33, hsRPB3; RPB31
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)