

Product datasheet for **GA103665**

Human mtRNA polymerase (POLRMT) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	POLRMT CRISPRa kit - CRISPR gene activation of human RNA polymerase mitochondrial
Symbol:	POLRMT
Components:	GA103665G1, mtRNA polymerase gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103665G2, mtRNA polymerase gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103665G3, mtRNA polymerase gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O00411</u>
Synonyms:	APOLMT; h-mtRPOL; MTRNAP; MTRPOL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)