

Product datasheet for **GA103669**

Human PON3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PON3 CRISPRa kit - CRISPR gene activation of human paraoxonase 3
Symbol:	PON3
UniProt ID:	<u>Q15166</u>
Synonyms:	paraoxanase-3; paraoxonase 3; serum paraoxonase/lactonase 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)