

Product datasheet for **GA103681**

Human POU4F3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | POU4F3 CRISPRa kit - CRISPR gene activation of human POU class 4 homeobox 3 |
| Symbol: | POU4F3 |
| UniProt ID: | Q15319 |
| Synonyms: | BRN3C; DFNA15 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)