

Product datasheet for **GA103688**

Human PPAR gamma (PPARG) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | PPARG CRISPRa kit - CRISPR gene activation of human peroxisome proliferator activated receptor gamma |
| Symbol: | PPARG |
| UniProt ID: | P37231 |
| Synonyms: | CIMT1; GLM1; NR1C3; PPARG1; PPARG2; PPARGgamma |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)