

Product datasheet for **GA103709**

Human PPP1R3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP1R3A CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 3A
Symbol:	PPP1R3A
UniProt ID:	Q16821
Synonyms:	GM; PP1G; PPP1R3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)