

Product datasheet for **GA103710**

Human PPP1R3C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP1R3C CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 3C
Symbol:	PPP1R3C
UniProt ID:	Q9UQK1
Synonyms:	PPP1R5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)