

Product datasheet for **GA103729**

Human Calcineurin A (PPP3CA) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP3CA CRISPRa kit - CRISPR gene activation of human protein phosphatase 3 catalytic subunit alpha
Symbol:	PPP3CA
UniProt ID:	<u>Q08209</u>
Synonyms:	CALN; CALNA; CALNA1; CCN1; CNA1; PPP2B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)