

## Product datasheet for **GA103734**

### Human PPP3R2 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PPP3R2 CRISPRa kit - CRISPR gene activation of human protein phosphatase 3 regulatory subunit B, beta
Symbol:	PPP3R2
UniProt ID:	<a href="#">Q96LZ3</a>
Synonyms:	PPP3RL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)