

Product datasheet for **GA103738**

Human Pancreatic Polypeptide (PPY) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPY CRISPRa kit - CRISPR gene activation of human pancreatic polypeptide
Symbol:	PPY
UniProt ID:	<u>P01298</u>
Synonyms:	PNP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)