

Product datasheet for **GA103739**

Human NPY4R activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NPY4R CRISPRa kit - CRISPR gene activation of human neuropeptide Y receptor Y4
Symbol:	NPY4R
UniProt ID:	P50391
Synonyms:	NPY4-R; PP1; PPYR1; Y4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)