

## Product datasheet for **GA103750**

### Human PRH1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | PRH1 CRISPRa kit - CRISPR gene activation of human proline rich protein HaellI subfamily 1                                                                                                                                                                                                              |
| Symbol:       | PRH1                                                                                                                                                                                                                                                                                                    |
| Components:   | GA103750G1, PRH1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA103750G2, PRH1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA103750G3, PRH1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>P02810</u>                                                                                                                                                                                                                                                                                           |
| Synonyms:     | Db-s; PA; PIF-S; Pr1/Pr2; PRH2                                                                                                                                                                                                                                                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)