

Product datasheet for **GA103764**

Human Protein Kinase A regulatory subunit I alpha (PRKAR1A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRKAR1A CRISPRa kit - CRISPR gene activation of human protein kinase cAMP-dependent type I regulatory subunit alpha
Symbol:	PRKAR1A
UniProt ID:	P10644
Synonyms:	ACRDYS1; ADOHR; CAR; CNC; CNC1; PKR1; PPNAD1; PRKAR1; TSE1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)