

Product datasheet for **GA103777**

Human PKC mu (PRKD1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRKD1 CRISPRa kit - CRISPR gene activation of human protein kinase D1
Symbol:	PRKD1
UniProt ID:	Q15139
Synonyms:	PKC-MU; PKCM; PKD; PRKCM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)