

Product datasheet for **GA103778**

Human PKC theta (PRKCQ) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRKCQ CRISPRa kit - CRISPR gene activation of human protein kinase C theta
Symbol:	PRKCQ
UniProt ID:	Q04759
Synonyms:	nPKC-theta; PRKCT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)