

Product datasheet for **GA103796**

Human MEK3 (MAP2K3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MAP2K3 CRISPRa kit - CRISPR gene activation of human mitogen-activated protein kinase kinase 3
Symbol:	MAP2K3
UniProt ID:	P46734
Synonyms:	MAPKK3; MEK3; MKK3; PRKMK3; SAPKK-2; SAPKK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)