

Product datasheet for **GA103810**

Human Protein C (PROC) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PROC CRISPRa kit - CRISPR gene activation of human protein C, inactivator of coagulation factors Va and VIIIa
Symbol:	PROC
UniProt ID:	P04070
Synonyms:	APC; PC; PROC1; THPH3; THPH4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)