

Product datasheet for **GA103835**

Human PR3 (PRTN3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRTN3 CRISPRa kit - CRISPR gene activation of human proteinase 3
Symbol:	PRTN3
UniProt ID:	P24158
Synonyms:	ACPA; AGP7; C-ANCA; CANCA; MBN; MBT; NP-4; NP4; P29; PR-3; PR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)