

Product datasheet for **GA103839**

Human Presenilin 1 (PSEN1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSEN1 CRISPRa kit - CRISPR gene activation of human presenilin 1
Symbol:	PSEN1
UniProt ID:	P49768
Synonyms:	AD3; FAD; PS-1; PS1; S182
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)