

Product datasheet for **GA103840**

Human Presenilin 2 (PSEN2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSEN2 CRISPRa kit - CRISPR gene activation of human presenilin 2
Symbol:	PSEN2
UniProt ID:	P49810
Synonyms:	AD3L; AD4; CMD1V; PS2; STM2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)