

## Product datasheet for **GA103843**

### Human PSG3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PSG3 CRISPRa kit - CRISPR gene activation of human pregnancy specific beta-1-glycoprotein 3
Symbol:	PSG3
UniProt ID:	<a href="#">Q16557</a>
Synonyms:	pregnancy specific beta-1-glycoprotein 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)