

Product datasheet for **GA103844**

Human PSG4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSG4 CRISPRa kit - CRISPR gene activation of human pregnancy specific beta-1-glycoprotein 4
Symbol:	PSG4
UniProt ID:	<u>Q00888</u>
Synonyms:	PSG9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)