

Product datasheet for **GA103850**

Human PSKH1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSKH1 CRISPRa kit - CRISPR gene activation of human protein serine kinase H1
Symbol:	PSKH1
UniProt ID:	<u>P11801</u>
Synonyms:	protein serine kinase H1; serine/threonine-protein kinase H1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)