

Product datasheet for **GA103851**

Human PSMA1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMA1 CRISPRa kit - CRISPR gene activation of human proteasome subunit alpha 1
Symbol:	PSMA1
UniProt ID:	P25786
Synonyms:	HC2; HEL-S-275; NU; PROS30
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)