

## Product datasheet for **GA103853**

### Human PSMA3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PSMA3 CRISPRa kit - CRISPR gene activation of human proteasome subunit alpha 3
Symbol:	PSMA3
UniProt ID:	<a href="#">P25788</a>
Synonyms:	HC8; PSC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)