

Product datasheet for **GA103854**

Human PSMA4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMA4 CRISPRa kit - CRISPR gene activation of human proteasome subunit alpha 4
Symbol:	PSMA4
UniProt ID:	P25789
Synonyms:	HC9; HsT17706; PSC9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)