

Product datasheet for **GA103856**

Human Proteasome subunit alpha type 6 (PSMA6) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMA6 CRISPRa kit - CRISPR gene activation of human proteasome subunit alpha 6
Symbol:	PSMA6
UniProt ID:	<u>P60900</u>
Synonyms:	IOTA; p27K; PROS27
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

