

Product datasheet for **GA103861**

Human Proteasome subunit beta type 4 (PSMB4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMB4 CRISPRa kit - CRISPR gene activation of human proteasome subunit beta 4
Symbol:	PSMB4
UniProt ID:	P28070
Synonyms:	HN3; HsN3; PROS-26; PROS26
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)