

## Product datasheet for **GA103868**

### Human PSMB10 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | PSMB10 CRISPRa kit - CRISPR gene activation of human proteasome subunit beta 10 |
| Symbol:       | PSMB10  |
| UniProt ID:   | <u>P40306</u>   |
| Synonyms:     | beta2i; LMP10; MECL1  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)         |

#### Product images:



[View online >](#)