

## Product datasheet for **GA103869**

### Human S4 (PSMC1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PSMC1 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, ATPase 1
Symbol:	PSMC1
UniProt ID:	<a href="#">P62191</a>
Synonyms:	P26S4; p56; S4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)