

## Product datasheet for **GA103872**

### Human Tbp7 (PSMC4) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | PSMC4 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, ATPase 4 |
| Symbol:       | PSMC4  |
| UniProt ID:   | <a href="#">P43686</a>   |
| Synonyms:     | MIP224; RPT3; S6; TBP-7; TBP7  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)              |

#### Product images:



[View online >](#)