

Product datasheet for **GA103874**

Human PSMC6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMC6 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, ATPase 6
Symbol:	PSMC6
UniProt ID:	P62333
Synonyms:	CADP44; HEL-S-73; p42; P44; SUG2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)