

Product datasheet for **GA103875**

Human PSMD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD1 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 1
Symbol:	PSMD1
UniProt ID:	Q99460
Synonyms:	P112; Rpn2; S1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)