

Product datasheet for **GA103877**

Human PSMD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD3 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 3
Symbol:	PSMD3
UniProt ID:	O43242
Synonyms:	P58; RPN3; S3; TSTA2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)