

Product datasheet for **GA103878**

Human PSMD4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD4 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 4
Symbol:	PSMD4
UniProt ID:	P55036
Synonyms:	AF; AF-1; ASF; MCB1; pUB-R5; Rpn10; S5A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)