

Product datasheet for **GA103879**

Human PSMD5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD5 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 5
Symbol:	PSMD5
UniProt ID:	Q16401
Synonyms:	S5B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)