

Product datasheet for **GA103880**

Human PSMD7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD7 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 7
Symbol:	PSMD7
UniProt ID:	P51665
Synonyms:	MOV34; P40; Rpn8; S12
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)