

## Product datasheet for **GA103881**

### Human PSMD8 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PSMD8 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 8
Symbol:	PSMD8
UniProt ID:	<a href="#">P48556</a>
Synonyms:	HEL-S-91n; HIP6; HYPF; Nin1p; p31; Rpn12; S14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)