

Product datasheet for **GA103882**

Human PSMD9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD9 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 9
Symbol:	PSMD9
UniProt ID:	O00233
Synonyms:	p27; Rpn4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)