

Product datasheet for **GA103883**

Human PSMD10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD10 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 10
Symbol:	PSMD10
UniProt ID:	O75832
Synonyms:	dj889N15.2; p28; p28(GANK)
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)