

Product datasheet for **GA103884**

Human PSMD11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD11 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 11
Symbol:	PSMD11
UniProt ID:	<u>O00231</u>
Synonyms:	p44.5; Rpn6; S9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)